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References

1. Hoxson S. Surface roughness parameters as predictors of exchange strength in bone: a critical analysis. *J Biomech* 2003;33(10):1297-303. (Ref. No. 7216)
2. Hoxson S, Hansson KN. The effect of limited lateral resolution in the measurement of implant surface roughness: a computer simulation. *J Biomed Mater Res A* 2005;75(2):472-7.
3. Hansson S, Norman M. The relation between surface roughness and interfacial shear strength for bone-anchored implants: A mathematical model. *J Biomech* 1999;32(8):829-36. (Ref. No. 7519)
4. Cooper LJ, Masuda T, Whitson SW, Yihelkissa F, Felton DA. Formation of autologous osteoblast colonies on machined, titanium oxide-grit-blasted, and plasma-sprayed titanium surfaces. *Int J Oral Maxillofac Implants* 1999;14(1):37-47.
5. Masuda K, Siva Lopez B, Holmberg J, Wronberg A, Arvinksson K. Attachment and proliferation of human oral fibroblasts to titanium surfaces blasted with TiO₂ particles: A scanning electron microscopic and histomorphometric analysis. *Clin Oral Implants Res* 1998;9(3):195-207.
6. Ericsson I, Johansson CB, Byström H, Norman MR. A histomorphometric evaluation of bone-to-implant contact on machine-prepared and roughened titanium dental implants: A pilot study in the dog. *Clin Oral Implants Res* 1994;5(4):285-6.
7. Grefröen K, Nims L, Hjerling-Hansen E, Jensen IS, Holman A. Histomorphometric and removal torque analysis for TiO₂-blasted titanium implants: An experimental study in dogs. *Clin Oral Implants Res* 1992;3(2):77-84.
8. Grefröen K, Wronberg A, Johansson C, Stenquist LT, Hjerling-Hansen E. Attachment of TiO₂-blasted, H₂O₂-etched, and machined implants: an experimental study with rabbits. *J Biomed Mater Res* 1995;29(1):123-31.
9. Rasmussen L, Kihlberg KE, Ten A. Effect of implant design and surface on bone regeneration and implant stability: an experimental study in the dog mandible. *Clin Impl Dent Rel Res* 2011;13(4). (Ref. No. 75411)
10. Kivallik CI, Holgers C, Widmark C, Senerby L, Wronberg A. Histological evaluation of the bone integration of TiO₂ blasted and turned titanium microimplants in humans. *Clin Oral Implants Res* 2011;12(2):128-34.
11. Abrahamson I, Berglundh T. Tissue characteristics at microblasted implants: an experimental study in dogs. *Clin Impl Dent Rel Res* 2006;8(3):167-8. (Ref. No. 7677)
12. Abrahamson I, Berglundh T, Moon ES, Lindhe J. Peri-implant tissues at submerged and non-submerged titanium implants. *J Clin Periodontol* 1999;26(7):600-7. (Ref. No. 7516)
13. Abrahamson I, Berglundh T, Selman S, Lindhe J. Tissue reactions to treatment shill: an experimental study in dogs. *Clin Impl Dent Rel Res* 2003;5(2):82-8. (Ref. No. 7622)
14. Alex A, Chang TL, Tekdemir I, Fousoo M. Biomechanical aspects of initial introsseous stability and implant design: a quantitative micro-morphometric analysis. *Clin Oral Implants Res* 2006;17(4):465-72.
15. Albury JP, Abrahamson I, Persson LG, Berglundh T. Spontaneous progression of peri-implantitis at different types of implants: An experimental study in dogs: A clinical and radiographic observation. *Clin Oral Implants Res* 2008;19(10):997-1002. (Ref. No. 79117)
16. Abrahamson I, Senerby L, Wronberg A. State of the art of oral implants. *Periodontol* 2000;20(8):473-86.
17. Al-Nawas G, Geil H. Three-dimensional topographic and histologic evaluation of dental implants by confocal laser scanning microscopy. *Clin Impl Dent Rel Res* 2003;5(3):176-83.
18. Arvidsson A, Sier BA, Wronberg A. The role of functional parameters for topographical characterization of bone-anchored implants. *Clin Impl Dent Rel Res* 2006;8(2):70-6.
19. Berglundh T, Abrahamson I, Albury JP, Lindhe J. Bone healing of implants with a fluoroide modified surface: an experimental study in dogs. *Clin Oral Implants Res* 2007;18(2):147-52. (Ref. No. 78678)
20. Berglundh T, Abrahamson I, Lindhe J. Bone reaction to biogranular (resorbable) and to implant: an experimental study in dogs. *J Clin Periodontol* 2005;32(9):925-32. (Ref. No. 78475)
21. Bettinelli D, Fermion LG, Lindhe J, Berglundh T. Bone tissue formation adjacent to implants placed in tooth extraction sockets: an experimental study in dogs. *Clin Oral Implants Res* 2006;17(4):351-58.
22. Boudrias VA, Sander-Pekerson S, Cantorighan JL, Brown C, Torres A, Nouri EE, et al. Immediate functional loading of single-tooth TiO₂ grit-blasted implant restorations: a controlled prospective study in a porcine model. Part I. Clinical outcome. *Clin Impl Dent Rel Res* 2007;9(4):197-206.
23. Boudrias VA, Wilhelmsson F, Jansen JA, Cunningham H, Brown C, Torres A, et al. Immediate functional loading of single-tooth TiO₂ grit-blasted implant restorations. A controlled prospective study in a porcine model. Part II. Histology and histomorphometry. *Clin Impl Dent Rel Res* 2007;9(4):207-46.
24. Carragolla D, Araujo M, Berglundh T, Albrektsson T, Lindhe J. Bone tissue reaction around implants placed in a compromised jaw. *J Clin Periodontol* 1999;26(3):429-35.
25. Carragolla D, Berglundh T, Araujo M, Albrektsson T, Lindhe J. Bone healing around implants placed in a jaw defect augmented with Bio-Ox. An experimental study in dogs. *J Clin Periodontol* 2003;27(1):819-885.
26. Cahn G, Gallat M, Nunez R, Borrajo JL, Alvarez JC, Vives L, et al. Scanning electron microscopic analysis of diode laser-treated titanium implant surfaces. *Plast Reconstr Surg* 2007;25(2):124-4.
27. Cooper LJ, Zhou Y, Takita J, Gao J, Abbin A, Holmberg J, et al. Fluoride modification effects on osteoblast behavior and bone formation at TiO₂ grit-blasted e.p. titanium endosseous implants. *Biomaterials* 2006;27(26):2826. (Ref. No. 7878)
28. Ellingsen JE, Johansson CB, Wronberg A, Holman A. Improved retention and bone-to-implant contact with fluoride-modified titanium implants. *Int J Oral Maxillofac Implants* 2001;16(5):459-66. (Ref. No. 78245)
29. Ellingsen JE, Lyngstad SA SF. Increasing biocompatibility by chemical modification of titanium surfaces. In: Ellingsen JE, Lyngstad SA SF, editors. *Bio-implant interface: Improving Biomaterials and Tissue Integration*. Boca Raton, Florida: CRC Press LLC; 2003. p. 323-40. (Ref. No. 7800)
30. Ericsson I, Lekholm U, Senerby L, Holman A. Soft tissue response to clinically contaminated and thereafter cleaned titanium surfaces: An experimental study in the rat. *Clin Oral Implants Res* 2004;15(4):370-3. (Ref. No. 75335)
31. Gao J, Fazioli RJ, Ambrose W, De Kolt H, Cooper LJ. The effect of hydrofluoric acid treatment of TiO₂ grit blasted titanium implants on adherent osteoblast gene expression in vitro and in vivo. *Int J Oral Maxillofac Implants* 2005;20(10):518-25.
32. Ise ZM, Schneider GG, Zaharias R, Senhold D, Stenford CM. Effects of fluoride-modified titanium surfaces on osteoblast proliferation and gene expression. *Int J Oral Maxillofac Implants* 2006;21(2):303-11. (Ref. No. 78731)
33. Isidor F. Loss of osseointegration caused by occlusal load of oral implants: A clinical and radiographic study in monkeys. *Clin Oral Implants Res* 1996;7(1):14-32.
34. Isidor F. Clinical printing and radiographic assessment in relation to the histologic bone level at oral implants in monkeys. *Clin Oral Implants Res* 1997;8(4):257-64.
35. Isidor F. Histological evaluation of peri-implant bone at implants subjected to occlusal overload or plaque accumulation. *Clin Oral Implants Res* 1997;8(1):1-9. (Ref. No. 7518)
36. Miskál C, Schneider GR, Zaharias R, Senhold D, Stenford C. Effects of implant surface microporosity on osteoblast gene expression. *Clin Oral Implants Res* 2005;16(6):652-6.
37. Mitzel L, Currie F, Jacobson M, Albrektsson T, Wronberg A. The effect of chemical and nanosporosity modifications on the early stages of osseointegration. *Int J Oral Maxillofac Implants* 2008;23(1):41-7.
38. Monje M, Lavella SF, Lyngstad SA SF, Rostad HL, Ellingsen JE. In vivo evaluation of osteogenic medium and bone mineral density at the surface of fluoride-modified titanium implants. *Biomaterials* 2008;29(26):3771-80.
39. Moon ES, Berglundh T, Abrahamson I, Lindhe J. The barrier between the keratinized mucosa and the dental implant: An experimental study in the dog. *J Clin Periodontol* 1999;26(10):658-63. (Ref. No. 75158)
40. Norton MR. An in vivo evaluation of the strength of an internal osseous interface compared to a bone joint interface in implant design. *Clin Oral Implants Res* 1997;8(4):293-8. (Ref. No. 75183)
41. Norton MR. In vivo evaluation of the strength of the cement implant-to-abutment joint in two commercially available implant systems. *J Prosthet Dent* 2000;83(5):607-71. (Ref. No. 72300)
42. O'Sullivan D, Senerby L, Mousdelli N. Microanalysis comparing the initial stability of five designs of dental implants in a human cadaver study. *Clin Impl Dent Rel Res* 2007;9(4):207-46.
43. Sponer P, Papakowits S, Hansson G, Brady K, Palmer RM, McDonald F. Cytokine release by osteoblast-like cells cultured on implant discs of varying alloy compositions. *Clin Oral Implants Res* 2002;13(6):623-30.
44. Thee A, Heng J, Zöllin G, Senerby L, Rasmussen L. Correlation of plaque growth factor release in vivo after placement of non-pure implants. accepted in *Clin Oral Implants Res* 2008;Work 4 in thesis by Weidner, 2008. ISBN 978-91-438-7982-4.
45. Zöllin G, Rodriguez VRL, Rizzo S, Bassetti C, Senerby L, Petroncini R, et al. Behavior of the bone-implant interface after push-in testing: a morphological study. *J Biomed Mater Res A* 2005;65(2):345-7.
46. Verstaeghe S, Walke JG, Næret I, Jensen JA. A histological evaluation of TiO₂ grit-blasted and Ca-P magnesium spinter coated implants placed into the trabecular bone of the goat. Part 1. *Clin Oral Implants Res* 2003;14(1):131-74.
47. Verstaeghe S, Walke JG, Næret I, Jensen JA. A mechanical evaluation of TiO₂ grit-blasted and Ca-P magnesium spinter coated implants placed into the trabecular bone of the goat. Part 1. *Clin Oral Implants Res* 2003;14(1):305-13.
48. Chaffee NR, Felton DA, Cooper LF, Palaniyandi U, Smith R. Prosthetic complications in an implant-retained mandibular overdenture population: initial analysis of a prospective study. *J Prosthet Dent* 2002;87(1):48-4.
49. Cooper LJ, Mottray JE, Guckes AD, Klee LB, Smith RG, Almgren C, et al. Five year prospective evaluation of mandibular overdentures retained by two microblasted, TiO₂ blasted remanufactured implants and retentive ball anchors. *Int J Oral Maxillofac Implants* 2008;23(4):696-704.
50. Cooper LF, Scarvia MS, Lang LA, Guckes AD, Mottray JE, Felton DA. Treatment of occlusal pain using Astra Tech implants and ball abutments to retain mandibular overdentures. *Int J Oral Maxillofac Implants* 1999;14(5):446-53. (Ref. No. 75155)
51. De Bruyn H, Besseler J, Rao F, Vanacker M. Clinical Outcome of Overdenture Treatment on Two Nitrogenized and Nonnitrogenated Astra Tech Microblasted maxillary Implants. *Clin Implant Dent Rel Res* 2008;10(4):310-11. (Ref. No. 8208:000103)
52. Kurita H, Sakai H, Uehara S, Kurahara K. Dental rehabilitation using an implant-carrying plate system in a severely resorbed edentulous maxilla: a case report. *Int J Oral Maxillofac Implants* 2005;20(2):111-20.
53. Makkonen TA, Holmberg S, Niemi L, Choon C, Tommila T, Felton J. A 5-year prospective clinical study of Astra Tech central implants supporting fixed bridges or overdentures in the edentulous mandible. *Clin Oral Implants Res* 1997;8(6):649-75. (Ref. No. 75181)

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